

# KPS2000 Series Piezo Jet Valve Upgraded Acoustics Semiconductor Optics Medical Stable And Consistent

### Basic Information

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity:
- Price: \$1-\$10000
- Packaging Details:
- Delivery Time: 5-60 days
- Payment Terms:



## **Product Specification**

- Type:
- KPS2000 Piezoelectric Valve Jetting System
- Autom

China

MingSeal

KPS2000

WOODEN

1

ISO CE

- Valve Dimension:
- Certification:

Automatic Grade:

- After-sales Service Provided:
- Warranty:
- Valve Weight:
- Highlight:

## Automatic

L/C,T/T,D/A,D/P,MoneyGram,Western Union

- 125.5\*75\*20mm(runner Assembled)
- CE Certificate, ISO, CE/ISO
  - Engineers Available To Service Machinery Overseas,Online Support,Video Technical Support,No Overseas Service Provided,Free Spare Parts
  - One Year

430g

Piezoelectric jet dispensing valve, Piezo jet dispensing valve, pick and place machine parts Valve



#### More Images



### KPS2000 Series Piezoelectric Jetting Valve (Upgraded)

KPS2000 bseries piezoelectric jetting valve is a new generation of high-speed, precision and non-contact jetting system which can deal with low, medium, high and ultra-high viscosity media; configurations such as UV type, corrosion-resistant type, PUR type and anti-corrosion type are available according to different application media.

#### Features and Advantages

Stable and Consistent Dispensing

The tappet guide seat is directly connected with the runner thread, and the coaxial and guiding design of the guide seat and the runner is added; the coaxiality of the tappet and the nozzle is improved to ensure continuously stable dispensing. The jet force of the valve body is increased by 20%, and

higher-viscosity glue (up to 500000cps) can be supported. Lower Usage Costs

With integrated design of nozzle seat, 1 sealing ring + 1 nozzle sealing seat + 1 installation tool are reduced to decrease maintenance and usage costs.

Faster Process Recovery Capability

Auxiliary enhanced version of ADJ calibration function and repeAfter routine cleaning and assembly, jetting can be done directly without ADJ calibration, and the process recovery time is less than 2 minutes.atability of 0.01mm.



#### Technical Specification

KPS2000 Plezoelectric Jetting Valve		KPC2000 Plezoelectrie Jeting Valve Controller	
Supply Pressure Range	0.1 8 Bar	Rated Power Supply	AC 220V±20%
Viscosity Range	Medium and high viscosity up to 500000mpas	Dispensing time and min.step 0.01ms-6000.00ms,0.01ms	
Maimum Jetting Freq	1200H2	Jetting Freq	1~1200HZ
Working Temp.Range	10 50°C	Display	TFT color screen
Corrosion Resistance	All water-based glue,organic	Parameters	10 sets
	solvents,weak acids and alkali	memory	
Valve Dimensions	125.5×40×20mm(w/o runner)	Reserved	Supporting STD.RS232/RS485 MODBUS
		Interface	RTU agreement
	125.5×75×20mm(runner	Controller	3.55KG
	assembled)	Weight	
Valve Weight	430g(with runner,w/o cable)	Controller	232×170×165mm
		Dimensions	
Supply Pressure Range	0.1 8 Bar	Rated Power	275W

Our Product Introduction

